

Reference values for serum Hepcidin-25

Reference ranges for serum hepcidin (nM) per 5-year age group for men and women in a reference population¹ as measured by weak cation exchange time-of-flight mass spectrometry (WCX-TOF MS)^{2,3,4}.

Reference levels for the WCX-TOF MS method are recalculated from those of our ELISA method¹, based on the regression line: (ELISA - 1.00)/1.52 = WCX-TOF MS that was derived from the results obtained by both methods for the same samples without hepcidin isoforms³. All values are determined using secondary reference material for hepcidin assays, which value is assigned by a primary reference material, allowing traceability to the internationally recognized Système International⁴. Results for heparine plasma, EDTA plasma, citrate plasma and serum do not differ from each other.

| Age, years | Men (n=1066) | | | | | Women (n=882)*, # | | | | |
|------------|-------------------|------|------------|---------------------|-------------|-------------------|------|------------|---------------------|-------------|
| | N | (%) | Median | 95% reference range | | N | (%) | Median | 95% reference range | |
| | | | | P2.5 | P97.5 | | | | P2.5 | P97.5 |
| 18-24 | 10 | (1) | 5.6 | 0.8 | 11.6 | 21 | (2) | 1.2 | < 0.5 | 6.6 |
| 25-29 | 16 | (2) | 5.2 | < 0.5 | 16.1 | 28 | (3) | 1.5 | < 0.5 | 7.0 |
| 30-34 | 18 | (2) | 4.4 | < 0.5 | 16.7 | 24 | (3) | 2.0 | < 0.5 | 13.8 |
| 35-39 | 22 | (2) | 3.8 | < 0.5 | 12.8 | 36 | (4) | 1.6 | < 0.5 | 10.4 |
| 40-44 | 19 | (2) | 6.4 | < 0.5 | 11.6 | 65 | (7) | 2.6 | < 0.5 | 16.0 |
| 45-49 | 76 | (7) | 3.6 | < 0.5 | 13.9 | 110 | (12) | 1.7 | < 0.5 | 9.4 |
| 50-54 | 106 | (10) | 4.2 | < 0.5 | 14.6 | 140 | (16) | 3.1 | < 0.5 | 15.1 |
| 55-59 | 173 | (16) | 4.6 | < 0.5 | 16.6 | 129 | (15) | 5.2 | < 0.5 | 14.3 |
| 60-64 | 179 | (17) | 4.7 | < 0.5 | 15.1 | 137 | (16) | 5.0 | < 0.5 | 18.2 |
| 65-69 | 186 | (17) | 5.6 | < 0.5 | 14.7 | 95 | (11) | 5.2 | < 0.5 | 15.0 |
| 70-74 | 133 | (12) | 5.2 | < 0.5 | 17.9 | 62 | (7) | 5.4 | < 0.5 | 25.5 |
| 75-79 | 99 | (9) | 4.0 | < 0.5 | 17.0 | 16 | (2) | 5.7 | 0.8 | 19.4 |
| 80-84 | 22 | (2) | 4.0 | 1.7 | 13.3 | 10 | (1) | 7.6 | < 0.5 | 12.7 |
| ≥85 | 7 | (1) | 7.2 | 1.7 | 13.5 | 9 | (1) | 3.9 | < 0.5 | 16.2 |
| All | 1066 (100) | | 4.7 | < 0.5 | 15.5 | 882 (100) | | 3.8 | < 0.5 | 15.4 |

*Pre-menopausal women (age <55; n=424)

Median = **2.1 nM**

P2.5 = 0.1 nM (<0.5 nM)

P97.5 = 13.0 nM

#Post-menopausal women (age ≥55; n=458)

Median = **5.2 nM**

P2.5 = 0.2 nM (<0.5 nM)


P97.5 = 16.5 nM

¹Galesloot TE, Vermeulen SH, Geurts-Moespot AJ, Klaver SM, Kroot JJ, van Tienoven D, Wetzel JF, Kiemeney LA, Sweep FC, den Heijer M, Swinkels DW. Serum hepcidin: reference ranges and biochemical correlates in the general population. *Blood* 2011; **117**: e218-25.

²Laarakkers CM, Wiegerink ET, Klaver S, Kolodziejczyk M, Gille H, Hohlbaum AM, Tjalsma H, Swinkels DW. Improved mass spectrometry assay for plasma hepcidin: detection and characterization of a novel hepcidin isoform. *PLoS ONE* 2013; **10**: e75518.

³Kroot JJ, Laarakkers CM, Geurts-Moespot AJ, Grebenchtchikov N, Pickkers P, van Ede AE, Peters HP, et al. Immunochemical and mass-spectrometry-based serum hepcidin assays for iron metabolism disorders. *Clin Chem* 2010; **56**: 1570-1579.

⁴Diepeveen LE et al. Provisional standardization of hepcidin assays: creating a traceability chain with a primary reference material, candidate reference method and a commutable secondary reference material. *Clin Chem Lab Med*. 2018 Nov **29**.



Reference values for [Hepcidin-25/Ferritin] and [TSAT/Hepcidin-25] ratios

Hepcidin values as given on page 1 should, like other hormones, be interpreted in the context of other indices of iron metabolism. For instance, in the absence of inflammation a low Hepcidin-25/ferritin ratio may be consistent with the presence of hereditary hemochromatosis or an iron loading anemia. On the other hand, a low transferrin saturation(TSAT)/Hepcidin-25 ratio may be consistent with Iron Refractory Iron Deficiency Anemia (IRIDA) due to a defect in the *TMPRSS6* gene.

WCX-TOF MS – [Hepcidin-25/ferritin] ratios

Men (n=1064)

Median = **28.2 pmol/µg**
 P2.5 = 3.1 pmol/µg
 P97.5 = 92.7 pmol/µg

Pre-menopausal women (age <55; n=424)

Median = **37.6 pmol/µg**
 P2.5 = 3.2 pmol/µg
 P97.5 = 176.4 pmol/µg

Post-menopausal women (age >=55; n=458)

Median = **42.7 pmol/µg**
 P2.5 = 9.6 pmol/µg
 P97.5 = 150.9 pmol/µg

WCX-TOF MS – [TSAT/Hepcidin-25] ratios

Men (n=1059)

Median = **6.9 %/nM**
 P2.5 = 1.6 %/nM
 P97.5 = 243.0 %/nM

Pre-menopausal women (age <55; n=422)

Median = **13.2 %/nM**
 P2.5 = 1.9 %/nM
 P97.5 = 312.9 %/nM

Post-menopausal women (age >=55; n=457)

Median = **5.4 %/nM**
 P2.5 = 1.4 %/nM
 P97.5 = 69.6 %/nM