
RBC Lysis

Research Use Only

- [Protocol A: RBC Lysis of Mouse Splenocytes](#)
- [Protocol B: RBC Lysis of Mouse Blood](#)
- [Protocol C: RBC Lysis of Human Peripheral Blood](#)
- [Protocol D: Lysis of Human Blood for Flow Cytometric Analysis](#)

Introduction

Prior to using lymphoid tissue cell suspensions for flow cytometric analysis and/or for *in vitro* functional assays, it is recommended to remove red blood cells. The eBioscience 1X RBC Lysis Buffer (Cat. No. [00-4333](#)) is formulated for optimal lysis of erythrocytes in single cell suspensions of mouse tissues such as spleen and human peripheral blood. The buffer contains ammonium chloride, which lyses red cells with minimal effect on lymphocytes. This step is not necessary when working with mouse thymus and lymph nodes.

When using human peripheral blood for flow cytometric analysis, a red cell lysing step is incorporated into the staining protocol.

For lysing human blood in bulk, a 10-20X volume of 1X RBC lysis buffer to total blood is recommended; however, an alternative method for removal of red blood cells when working with human blood is to use isolation of PBMC by Ficoll-Hypaque™.

Protocol A: RBC Lysis of Mouse Splenocytes

Materials

- 1X PBS:
 - 80.0g NaCl 80.0g NaCl
 - 11.6g Na₂HPO₄
 - 2.0g KH₂PO₄
 - 2.0g KCl
 - DI H₂O up to 10.0 L
 - pH to 7.0
- eBioscience 1X RBC Lysis Buffer (Cat. No. [00-4333](#))
- 50ml conical tubes

Instruments Pipettes and pipettors

- Centrifuge
- Hemacytometer and microscope

Experiment Duration

- 20 minutes

RBC Lysis

Research Use Only

Experimental Procedure

1. Harvest mouse spleen and prepare a single cell suspension.
2. Pellet the cells by centrifugation (300-400xg) at 4°C and aspirate the supernatant.
3. Resuspend the pellet in 5ml/spleen of Lysis Buffer.
4. Incubate at room temperature for 4-5 minutes with occasional shaking (we have performed this step on ice successfully too).
5. Stop the reaction by diluting the Lysis Buffer with 20-30ml of 1X PBS.
6. Spin the cells (300-400xg) at 4°C and resuspend the pellet in the appropriate buffer for use in the next step of your experimental procedure.
7. Perform a cell count at this time.

Note: In general, a small number of residual red cells does not interfere with proliferation assays or can be gated out from flow cytometric analysis. However, if required, a second round of lysis can be performed.

Protocol B: RBC Lysis of Mouse Blood

Materials

- 1X PBS:
 - 80.0g NaCl 80.0g NaCl
 - 11.6g Na₂HPO₄
 - 2.0g KH₂PO₄
 - 2.0g KCl
 - DI H₂O up to 10.0 L
 - pH to 7.0
- eBioscience 1X RBC Lysis Buffer (Cat. No. [00-4333](#))
- 50ml conical tubes

InstrumentsPipettes and pipettors

- Centrifuge
- Hemacytometer and microscope

Experiment Duration

- 20 minutes

RBC Lysis

Research Use Only

Experimental Procedure

1. Add 10 ml of RBC lysis buffer per 1 ml of mouse blood.
2. Incubate at room temperature for 4-5 minutes with occasional shaking (we have performed this step on ice successfully too).
3. Stop the reaction by diluting the Lysis Buffer with 20-30 ml of 1X PBS.
4. Spin the cells (300-400xg) at 4°C and resuspend the pellet in the appropriate buffer for use in the next step of your experimental procedure.
5. Perform a cell count at this time.

Note: In general a small number of residual red cells does not interfere with the proliferation assays and can be gated out from flow cytometric analysis. However, if required, a second round of lysis can be performed.

Protocol C: RBC Lysis of Human Peripheral Blood

Materials

- 1X PBS:
 - 80.0g NaCl 80.0g NaCl
 - 11.6g Na₂HPO₄
 - 2.0g KH₂PO₄
 - 2.0g KCl
 - DI H₂O up to 10.0 L
 - pH to 7.0
- eBioscience 1X RBC Lysis Buffer (Cat. No. [00-4333](#))
- 50ml conical tubes

InstrumentsPipettes and pipettors

- Centrifuge
- Hemacytometer and microscope

Experiment Duration

- 20 minutes
1. Add 10 ml of lysis buffer per 1ml of human blood.
 2. Incubate for 10 minutes at room temperature (no more than 15 minutes).
 3. Stop the reaction by diluting the Lysis Buffer with 20-30 ml of 1X PBS.
 4. Spin the cells (300-400xg) at 4°C and resuspend the pellet in the appropriate buffer for use in the next step of your experimental procedure.
 5. Perform a cell count at this time.

Note: In general a small number of residual red cells does not interfere with the proliferation

RBC Lysis

Research Use Only

assays and can be gated out from flow cytometric analysis. However, if required, a second round of lysis can be performed.

Protocol D: Lysis of Human Blood for flow cytometric analysis

When using human peripheral blood for flow cytometric analysis, the necessary red cell lysing step is incorporated into the staining protocol. The staining protocol may be found at:
<http://www.ebioscience.com/ebioscience/appls/FCS.htm>