

In situ hybridization of DIG-labelled RNA probes
Prepared by Carin Lundmark

Chemicals

Acetic anhydride	SIGMA, A 6404
Albumin, bovine	SIGMA, A-9647
AquaPerm Mounting medium, 290 ml	Immunon, 484990
5-Bromo-4-chloro-3-inolyl-phosphate (BCIP)	B.M., 1585002
Blocking Reagent	B.M., 1 096 176
Dextran sulphate	SIGMA, D 8906
Diethyl pyrocarbonate (DEPC)	SIGMA, D 5758
Dimethylformamide, 100%	BDH, Cat.No. 10322
Ethanol, 95% & 99.5%	Kemetyl
Ficoll	SIGMA, F-2637
Formamide	MERCK, 1.09684
Glycine pa.	MERCK, 1.04201
Hydrochloric acid fuming	MERCK, 1.00317
4-Nitro blue tetrasodium chloride (NBT)	B.M, 1087479
Magnesium chloride hexahydrat pa.	MERCK, 1.05833
Maleic acid for syntesis	MERCK Sch., 800380
Mixed Bed Resin AG, 501-X8 Resin 20-50 mesh	BIORAD
Paraformaldehyde for syntesis	MERCK Sch., 818715
Pertex mounting medium	HISTOLAB Prod., 00811
Potassium chloride pa. (KCl)	MERCK, 4936
Potassium dihydrogen phosphate pa. (KH ₂ PO ₄)	MERCK, 4873

Proteinase K	SIGMA, P-2308
Polyvinyl pyrrolidone (PVP)	SIGMA, P-5288
Ribonuclease A (RNase A)	SIGMA, R-6513
Rabbit anti-DIG-AP, Fab fragments, rabbit F(ab')	DAKO, D5105
Ribonucleic acid transfer, from Bakers Yeast, Type X	SIGMA, R-9001
Sodium chloride pa. (NaCl)	MERCK, 1.06404
Sodium hydrogen phosphate, anhydrous pa. (Na ₂ HPO ₄)	MERCK, 1.06586
Titriplex III pa. EDTA	MERCK, 1.08418
Triethanolamine pa.	MERCK, 1.08379
Tri-sodiumcitrate dihydrate	MERCK, 1.06448
Tris, Molecular Biology Certified	IBI, IB70142
Triton X-100	SIGMA, X-100
Xylene pa.	MERCK, 1.08681
Digoxigenin(DIG)-labelled mRNA probe (riboprobe), prepared according to Boehringer Mannheim protocol.	

Stock and working solutions

50X Denhardt's Reagent

		<i>Final conc.</i>
500 mg	Ficoll	1%
500 mg	Polyvinylpyrrolidone (PVP)	1%
500 mg	RNase-free bovine serum albumin	10 mg/ml

Add DEPC-H₂O to final volume of 50 ml.

Hybridization buffer

		<i>Final conc.</i>
10 ml	Formamide, deionized	50%
2 ml	3 M NaCl	0.3 M

0.2 ml	2 M Tris-HCl, pH 8	10 mM
0.04 ml	0.5 M EDTA, pH 8	1 mM
1 ml	10 mg/ml tRNA	0.5 mg/ml
4 ml	0.5 g/ml DS	100 mg/ml
0.4 ml	50X Denhardt's reagent	1X
2.36 ml	DEPC-H ₂ O	

Aliquot in eppendorf tubes. Hybridization buffer can be stored at -20°C for several months.

4% Paraformaldehyde, pH 7.0; 200 ml

8 g	Paraformaldehyde
600 µl	5 M NaOH (in DEPC-H ₂ O)
20 ml	10X PBS (in DEPC-H ₂ O)

Triethanolamine-HAc; 200ml

2.66 ml	Triethanolamine
200 ml	DEPC-H ₂ O

Submerge slides, add 500 ml acetic anhydride. Mix the two solutions by moving slides in the solution.

Prehybridization mixture; 200ml

20 ml	20X SSC (in DEPC-H ₂ O)
100 ml	Formamide, deionized
80 ml	DEPC-H ₂ O

Hybridization solution

Add riboprobe to hybridization buffer: 0.5 ng probe/µl . Place the solution for 3 min on heating block 92-95°C.

Cool immediately on ice.

100 µl hybridization solution/glass is applied per glass

Hybridization solution is prepared freshly just before the hybridization step.

Buffer 1; 800 ml

		<i>Final conc.</i>
80 ml	1 M Tris-HCl, pH 7,5	100 mM
40 ml	3 M NaCl	150 mM

Add autoclaved, deionized water to final volume of 800 ml.

Buffer 1 with Triton; 800 ml

10 ml Buffer 1
800 µl Triton X-100

Add Triton X-100 on the inner wall of the test tube and dilute with Buffer 1.
Leave solution on a shaking platform for a few minutes.

RNase buffer; 600 ml

100 ml	3 M NaCl	<i>Final conc.</i> 500 mM
6 ml	1 M Tris-HCl, pH 8	10 mM
1.2 ml	0.5 M EDTA, pH 8	1 mM

Add autoclaved, deionized water to final volume of 600 ml.

20 mg/ml RNase; 150 ml

3 mg RNase A
150 ml RNase buffer

Antibody solution; 3 ml

Dilute the Anti-DIG-AP antibody 1:1000 in Buffer 1 with Triton. Apply 100 µl/glass.

Buffer 2; 200 ml

		<i>Final conc.</i>
20 ml	1 M Tris-HCl, pH 7.5	100 mM
10 ml	3 M NaCl	150 mM
200 ml	Triton X-100	0.1 %
20 ml	10% Blocking Reagent	1 %

Add autoclaved, deionized water to final volume of 200 ml.

Buffer 3; 200 ml

		<i>Final conc.</i>
20 ml	1M Tris-HCl, pH, 9.5	100 mM
5 ml	4 M NaCl	100 mM
10 ml	1 M MgCl ₂	50 mM

Add autoclaved, deionized water to final volume of 200 ml.

Substrate solution; 3 ml

3 ml	Buffer 3
13.5 µl	NBT
10.5 µl	BCIP

TE 8; 200 ml

		<i>Final conc.</i>
1 ml	2 M Tris-HCl, pH 8	10 mM
400 ml	0.5 M EDTA, pH 8	1 mM

Add autoclaved, deionized water to 200 ml.

Paraffin sections protocol:

Day 1

1. Xylene	2 x 5 min	
-----------	-----------	--

2. abs EtOH	Rinse	
3. abs EtOH	5 min	
4. 95% EtOH	Rinse	
5. 70% EtOH	Rinse	
6. 1X PBS	Rinse	
7. 0.2M HCl	15 min	
8. 1X PBS	2 x 5 min	
9. 20 µg/ml Proteinase K	30 min, 37°C	Use humid box.
10. 0.2% Glycine	Rinse	
11. 4% Paraformaldehyde	5 min	
12. Triethanolamine. When slides are in the dish: Add 500 µl acetic anhydride. Mix solutions by dipping slides several times.	10 min	Shaking platform
13. 1X PBS	Rinse	
14. 1X PBS	5 min	
15. Prehybridization mixture	60 min, 42°C	Water bath
16. Hybridization solution	Over night, 42°C	Draw a fat-ring around the section. Apply 100 µl/glass. Use humid box.

Day 2

17. 4X SSC	Rinse	Rinse each probe separately
18. 2X SSC	2 x15 min, 37°C	Oven with shake board
19. 1X SSC	2 x15 min, 37°C	Oven with shake board
20. 20 µg/ml RNase	30 min, 37°C	Special oven
21. RNase buffer	Rinse, 37°C	
22. RNase buffer	5 min, 37°C	Special oven
23. 2X SSC, 50% formamide	2 x 30 min, 52°C	Oven with shake board
24. Buffer 1	Rinse	
25. Buffer 1	5 min	
26. Buffer 2	30 min	Shaking platform.
27. Antibody solution	30 min, 37°C	Oven. Apply 100 µl/glass.
28. Buffer 1	Rinse	

29. Buffer 1	5 min	
30. Buffer 3	2 min	
31. Substrate solution	15 min to 3 hours (in the dark).until desired colour has eveloped. Examine sample by microscope.	Apply 100 µl/glass.
32. TE8	Rinse	

33. Mount sections with AquaPerm. Apply 1-3 drops of AquaPerm over the section. Leave section to dry: 1 hour at room temperature or 20-30 minutes at 40-60°C in oven. Make sure that the surface is completely dry and hard before you mount and coverslip sections with Pertex.

Cryo-sections protocol:

Day 1

1. 95% EtOH	5 min	
2. 70% EtOH	Rinse	
3. 1X PBS	2 x 5 min	
4. 20 µg/ml Proteinase K	30 min, room temp	Humid box
10. 0.2% Glycine	Rinse	
11. 4% Paraformaldehyde	5 min	
12. Triethanolamine. When slides are in the dish: Add 500 µl acetic anhydride. Mix solutions by dipping slides several times	10 min	Shaking platform
13. 1X PBS	Rinse	
14. 1X PBS	5 min	
15. Prehybridization mixture	60 min, 42°C	Water bath
		Draw a fat-ring around the section
16. Hybridization solution	Over night, 42°C	Apply 100 µl/glass. Use humid box.

Day 2

17. 4X SSC	Rinse	Rinse each probe separately.
18. 2X SSC	2 x 15 min, 37°C	Oven with shake board
19. 1X SSC	2 x 15 min, 37°C	Oven with shake board

20. 20 µg/ml RNase	30 min, 37°C	Special oven
21. RNase buffer	Rinse, 37°C	
22. RNase buffer	5 min, 37°C	Special oven
23. 2X SSC, 50% formamide	2 x 30 min, 52°C	Oven with shake board
24. Buffer 1	Rinse	
25. Buffer 1	5 min	
26. Buffer 2	30 min	Shaking platform
27. Antibody solution	30 min, 37°C	Oven. Apply 100 µl/glass.
28. Buffer 1	Rinse	
29. Buffer 1	5 min	
30. Buffer 3	2 min	
31. Substrate solution	15 min to 3 hours (in the dark) until desired colour has developed. Examine sample by microscope.	Apply 100 µl/glass.
32. TE8	Rinse	

33. Mount sections with AquaPerm. Apply 1-3 drops of AquaPerm over the section. Leave section to dry: 1 hour at room temperature or 20-30 minutes at 40-60°C in oven. Make sure that the surface is completely dry and hard before you mount and coverslip sections with Pertex.

Tissue preparations

Fixation

Keep the time between killing of animal and fixation as short as possible, for best preservation of morphology.

Place tissue in freshly prepared paraformaldehyde, 4%.

Fix the tissue over night (not more than 24 hours) at 4°C.

Decalcification

Discard paraformaldehyde solution. (Rinsing tissue is not necessary.)

Cover the tissue with EDTA-solution (~12.5%, pH 7) for decalcification.

Keep at 4°C with gentle shaking.

Change EDTA-solution every 3-4 days.

Check end-point with x-ray.

After decalcification: Rinse the samples in 1X PBS, some hours to over night, on shaking platform.

Dehydration and paraffin infiltration

Use Vacuum Infiltration Processor (V.I.P.). Program 1 is suitable fo rat jaws.

Keep blocks in refrigerator until use. Long storage times up to years are possible.

Freezing tissue

After completed decalcification , leave tissue in 10% sucrose over night. Freeze the tissue in liquid nitrogen or hexane on dry ice. Store in -70°C freezer.

Sectioning

Paraffin

Wear gloves during sectioning procedures.

Cut 7 μm thick sections from the paraffin-embedded material.

The water bath should contain DEPC-H₂O.

Place sections on silan-coated microscope slides.

Leave sections to dry at 56°C in oven over night.

Keep section in refrigerator until use; for maximally 1 month.

Cryo

Mount tissue on a preholder, use cryomount to stabilize the tissue.

Cut sections at 7 μm . Mount section on silan coated slides.

Leave section to dry in 56°C over night. Store sections at 4°C until use.